
System Installation

l'installation du système / Systemvorbereitung / Configuración del sistema /
 Configuração do sistema / Sistema installazione / Opstelling van het systeem /
 Opsætning af systemet / Systemoppsett / Asennus järjestelmä / Systeminstallationen /
 Sistem Kurulumu

VERSAJET[®] II EXACT VERSAJET[®] II PLUS

Hydrosurgery System Handpieces

Do not use contents if package is opened or damaged.
Single use Handpiece Assembly — Do not re-sterilize.

Indications for Use

The VERSAJET II Hydrosurgery System is intended for wound debridement (acute and chronic wounds, burns), soft tissue debridement, and cleansing of the surgical site in applications that, in the physician's judgment, require sharp debridement and/or pulsed lavage irrigation.

Device Description

Producing a high-pressure stream of sterile saline, the VERSAJET II Hydrosurgery System creates a localized vacuum in the operating window of the handpiece.

As the handpiece tip is applied parallel/tangential to the wound surface, various tissue types can be targeted for debridement. The design of the handpiece allows precise removal of tissue on large areas and irregular, contoured surfaces such as web spaces and facial structures.

Warnings

- This device should be used with particular care in patients with hemophilia or other blood clotting disorders and in patients receiving anti-coagulant medication.
- VERSAJET II EXACT and PLUS handpieces are designed for connection only to the VERSAJET II Hydrosurgery System power console (part number (p/n) 66800039).
- Increasing settings on the device console will lead to more aggressive tissue removal. Use caution near sensitive tissue, such as neurovascular bundles.
- This device can cut soft tissue. Apply only to tissue intended to be cut or to remove debris from the wound.
- VERSAJET II is intended primarily for use in the operating room environment. Only 45° VERSAJET II EXACT handpieces (p/n 66800041 and 66800042) are suitable for use outside the operating room. Attention to universal infection control procedures should be applied when using the device outside the controlled environment of the operating room.
- If VERSAJET II EXACT 45° handpieces are used outside the operating room, ensure the floor immediately surrounding the treatment area is covered and any splashes cleaned after treatment is complete.
- The 15°/14mm VERSAJET II EXACT handpiece (p/n 66800040) and all VERSAJET II PLUS handpieces (p/n 66800043, 66800044 and 66800045) should not be used outside the operating room due to potential excessive misting and/or spraying.
- VERSAJET II PLUS handpieces provide more aggressive debridement and excision when compared to VERSAJET II EXACT handpieces. VERSAJET II PLUS is appropriate for wounds that require aggressive, yet selective, removal of tissue that are tough, heavily necrotic, contaminated or burnt. Users should be aware that just as with any sharp instrument, care and attention must be maintained while using VERSAJET II PLUS handpieces near fragile or sensitive structures.

Precautions

- Always begin debridement procedures at lower power settings and increase as necessary to the optimal power setting for the type of tissue being debrided and to avoid unintended tissue removal.
- In order to avoid unwanted procedural delays, ensure the system is fully operational prior to administration of anesthesia.
- When used on wounds where bone, tendon or other dense tissue may be encountered during the debridement procedure, excessive spraying and/or misting may occur due to the interruption of the stream of sterile saline by hard tissue.
- Spraying and/or misting is more frequent at lower power settings due to less pressure. Spraying and/or misting may be reduced by keeping the waste evacuation tube straight.
- As with all surgical procedures, the VERSAJET II operator and clinical personnel should follow universal precautions for infection control (including the use of surgical gloves, facemask that covers the mouth and nose, protective eye goggles, protective clothing, and anti-slip shoe covers).
- For optimal results when debriding hard or leathery eschar resulting from burn injuries, it is recommended to first debride the eschar using sharp debridement techniques followed by the use of VERSAJET II to complete debridement or excision of the wound.
- It is recommended that prior to clinical use of VERSAJET II devices and handpieces, all operators should be trained in the proper use.

Smith & Nephew has developed a training program. Contact your local market representative for details concerning VERSAJET II training.

- Do not add antibiotics, medication or other fluids to the saline irrigation source. Additives may change fluid dynamics and negatively affect product performance.
- Do not allow saline bag to empty; this could allow air to enter the supply tube. Air in the supply tube will temporarily lower device efficiency and may require re-priming of the system.
- Connecting the waste evacuation tubing hose (or any container connected to it) to a vacuum source is not recommended and may increase the aggressiveness of tissue removal.
- Monitor fluid level of waste container and empty as needed.
- Subsequent debridement procedures may be necessary for complex or highly contaminated wounds.
- Do not touch the high pressure jet in the operating window of the handpiece.
- Use only sterile saline solution with this device.
- Examine all components before use. If you believe a component to be faulty, damaged, or suspect, DO NOT USE. Contact your local Smith & Nephew VERSAJET II representative.
- Each VERSAJET II disposable handpiece is intended for SINGLE USE ONLY. DO NOT RE-STERILIZE. Discard after use. VERSAJET II handpieces are not designed to withstand the rigors of reprocessing or re-sterilization; device performance will be compromised and sterility can not be assured. Refer to our company position in the VERSAJET II Hydrosurgery System User Manual or visit our web site: www.versajet.info
- Pre-heating saline prior to use with VERSAJET II is not recommended. Due to the use of high pressure, some heating of the saline will occur during use.

Special Storage and Handling Instructions

Do not store product above 125°F (52°C) or below -40°F (-40°C). For further information, consult the VERSAJET II Hydrosurgery System User Manual.

Instructions for Use

1. Remove the pouch from carton. Inspect pouch to ensure seals are intact and pouch is undamaged. Open the pouch ensuring that the sterility of the inner pouch is not compromised. Note: Chart labels are on the outer pouch.
2. Inspect inner pouch to ensure seals are intact and pouch is not damaged. Open inner pouch, remove the sterile contents and place securely in the sterile field. Avoid tangling and knotting of hoses.
3. Remove white handpiece from plastic tray and place in sterile field. Do not remove pump cartridge (orange handle) from clam shell tray.
4. Remove the white paper tape from coiled tubing. Uncoil the various tubing lines. Maintain aseptic technique for the white handle portion of the handset and several feet of tubing to allow access to the surgical site. Hand off the saline spike line, orange pump cartridge assembly in clam shell tray and clear evacuation line to the Circulator for final system set up.
5. The Circulator should remove the orange pump cartridge from the clam shell and insert into the handpiece connection port located on the front of the console until fully seated and then turn clockwise to the 3 o'clock position. When locked correctly, the circular light surrounding the connection port should illuminate green.
6. CAUTION: Only insert the saline bag spike AFTER the orange pump cartridge has been securely locked in the console. Failure to do so may result in fluid leakage from the pump cartridge.
7. Remove sterile cover from bag spike and insert into sterile saline supply bag. Ensure that the clip on the saline inflow line is open.
8. Attach end of waste evacuation tubing (blue tip) to waste container. DO NOT connect to a port containing a filter or to the port labeled "VACUUM". You must ensure that there is an additional open port on the waste containers lid. Ensure there are no kinks or other external obstructions in saline supply, high-pressure and waste evacuation hoses.
9. Insert the foot pedal cable into the connector on the front of the console until the cable locks in place. The red dot on the plug should be aligned with the red dot on the console. Position the foot pedal for convenient access.
10. Insert the female three-prong plug of the power cord into the back of the console and the other end into an electrical outlet. Press the main power switch located on the front of the console.
11. Remove the protective cover from the tip of handpiece.
12. While keeping the handpiece at a safe distance, set the console power setting to 10. Press the foot pedal and keep the system running until saline reaches the distal tip of the handpiece. An audible hissing sound and a visible saline jet aimed down the evacuation line, indicates the system is primed. Priming takes approximately 30 seconds. Release the foot pedal and reduce the power setting to 1 before use.

Rx only

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician

STERILE EO

Sterilized using ethylene oxide

Do not use if package is damaged

2 Do not reuse

LATEX

This product does not contain natural rubber latex

13. CAUTION: Once the system has been primed with saline, do not allow the saline bag to empty. An empty bag can cause air to enter the system and require re-priming of the system. Tubing should be clamped when changing saline bags.

14. Begin debridement procedures starting at the lowest power setting and increase as necessary to the optimal power setting for the type of tissue being excised or debrided.

15. If the device tip becomes blocked with foreign matter, a reduction in device efficiency or the presence of spray from the tip may result. To eliminate the obstruction, remove the handpiece from the wound site, release the foot pedal and remove the obstruction with forceps. Do not touch the opening in the high-pressure jet with forceps. Once the obstruction is removed, press the foot pedal and check for steady stream of sterile saline flow. If the obstruction is not completely removed, repeat procedure or check that the evacuation line is not pinched by forceps or being stepped on, or that the collection container is full.

16. After completing the procedure, turn off the console main power switch. Disconnect the handpiece assembly from the power console by turning the orange pump assembly counterclockwise to the 12 o'clock position. Remove the pump cartridge by pulling it straight out. Discard handpiece assembly in accordance with your healthcare facility's standard guidelines for biohazardous waste disposal.

17. Follow your healthcare facility standard procedures for decontaminating surgical equipment. The following are the manufacturer's recommendations for console decontamination: Wipe all surfaces of the console with paper towels or a cloth soaked in a disinfectant solution. Sporidicid® Disinfectant (1.12% glutaraldehyde 1.93% phenol/phenate) or an equivalent high-level disinfectant is recommended. Dilute the disinfectant solution according to the manufacturer's instructions. After disconnecting the foot pedal and power cable from the power source, wipe down all exposed surfaces of both components in accordance with the guidance for the console. Dispose of towels, gloves and gown in accordance with your healthcare facility's standard guidelines for biohazardous waste disposal.

Product Availability

Part Number	Size
66800040	VERSAJET II EXACT 15°/14mm
66800041	VERSAJET II EXACT 45°/14mm
66800042	VERSAJET II EXACT 45°/8mm
66800043	VERSAJET II PLUS 15°/14mm
66800044	VERSAJET II PLUS 45°/14mm
66800045	VERSAJET II PLUS 45°/8mm

Operating Characteristics

Values provided in the table are typical values and are provided as guidance only. Actual values may vary.

Power Setting	Fluidjet Flow Rate mL/min	VERSAJET II EXACT Fluid Pressure psi/bar	VERSAJET II EXACT Fluidjet Velocity mph/kph	VERSAJET II PLUS Fluid Pressure psi/bar	VERSAJET II PLUS Fluidjet Velocity mph/kph
1	90	1787/123	333/536	2100/145	414/666
2	109	2789/192	396/637	3100/214	501/806
3	123	3799/262	447/719	4300/297	567/912
4	138	4816/332	501/807	5500/379	632/1018
5	160	5924/408	590/949	6700/462	736/1184
6	166	7080/488	603/971	7900/545	764/1230
7	181	8180/564	658/1058	9200/634	830/1336
8	195	9297/641	709/1140	10400/717	896/1443
9	209	10421/719	759/1222	11700/807	962/1549
10	220	11535/795	810/1304	12800/883	1012/1628

Technical Assistance

Contact your local market representative for technical assistance.

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
China

Smith & Nephew Ltd.
03-06,10 Floor,
Headquarters Building,
168 Xizang Road (M),
HuangPu District, Shanghai, China
200001
Tel: +86-021-63503100
Fax: +86-021-63501219
info.china@smith-nephew.com

Denmark

Smith & Nephew A/S
Slotsmarken 14
2970 Hoersholm
Danmark
Tel: +45 4580 6100
Fax: +45 4580 6151
sn.dk@smith-nephew.com

Finland

Smith & Nephew Oy
Äyritie 12 C
01510 Vantaa
Finland
Tel: +358 (0)207 866 300
Fax: +358 (0)207 866 320
Tilaukset: +358 (0)207 866 333
etunimi.sukunimi@smith-nephew.com

France

Smith & Nephew SAS
Espace Novaxis
25 Blvd. Marie et Alexandre Oyon
72019 Le Mans Cedex 2
France
Tel: (33) 02 43 83 23 82
(33) 02 43 83 24 14
Fax: (33) 02 43 83 23 83
scw@smith-nephew.com

Germany

Smith & Nephew GmbH
Osterbrooksweg 71
D-22869 Schenefeld
Tel: +49 (040) 87 97 44 -0
Fax: +49 (040) 87 97 44-375
info@smith-nephew.com

Hong Kong

Smith & Nephew Ltd
Unit 813-816, 8th Floor
Delta House
3 On Yiu Street
Shatin NT
Hong Kong
Tel: +852 2648 7700
Fax: +852 2648 7282
hongkong@smith-nephew.com

India

Smith & Nephew Healthcare Pvt. Ltd.
B-203, Dynasty Business Park,
J. B. Nagar, Andheri Kurla Road,
Andheri East,
Mumbai - 400 059
India
Tel: +91 22 40055090
Fax: +91 22 40055089
india.admin@smith-nephew.com

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Singapore

Smith & Nephew Pte Ltd
2 Jurong East Street 21
IMM Building #04-32B
Singapore 609601
Tel: 65-62700552
Fax: 65-62726698
singapore.office2@smith-nephew.com

South Africa

Smith & Nephew (Pty) Ltd.
PO Box 92
Pinetown 3600
South Africa
Tel: +27 (0)31 242 8111
Fax: +27 (0)31 242 8120
lorraine.speed@smith-nephew.com

South Korea

Smith & Nephew Ltd.
ASEM Tower(13F)
Samsung-dong
Gangnam-gu
Seoul
(Seoul 135-798 Korea)
Tel: (02)6001-7575
Fax: (02)6001-7560

Spain

Smith & Nephew S.A.
Fructuós Gelabert, 2 y 4
08970 Sant Joan Despí
Barcelona
España
Tel: +34 93 373 7301
Fax: +34 93 373 7453
maria.moreno@smith-nephew.com

Sweden

Smith & Nephew AB
Kråketorpsgatan 20
Box 143
431 22 Mölndal
Sverige
Tel: 031-746 58 00
Fax: 031-87 05 32

Switzerland

Smith & Nephew AG
Glutz Blotzheim-Strasse 1
4502 Solothurn
Switzerland
Tel: +41 32 624 56 60
Fax: +41 32 624 56 81
swiss.info@smith-nephew.com

Taiwan

Smith & Nephew Taiwan
9F-3, No. 130, Sec. 2
Jhongsiao East Road
Taipei City, Taiwan (R.O.C.)
Tel: 886 2 3322 5020
Fax: 886 2 3322 5021

Thailand

Smith & Nephew Ltd
344/3 Soi Rongrien Yepun
Rama IX Road
Kwang Bangapi
Khet Huay Kwang
Bangkok 10320
Thailand
Tel: +66 2 719 6221-27
Fax: +66 2 719 6220
snlthai@loxinfo.co.th

Turkey

Baral Tibbi Malzemeler Ltd Sti
Nilgun Sok. No. 12/12 06680
Kavaklidere
Ankara
Turkey
Tel: +90 312 427 1167
bahadiralev@baral.com.tr

United Arab Emirates

Smith & Nephew FZE
Level 4,
National Bank of Abu Dhabi Bldg.
Dubai Health Care City
P.O.Box 9715
Dubai
United Arab Emirates
Tel: +971 4 4299111
Fax: +971 4 4299100
snfze@emirates.net.ae

United Kingdom

Smith & Nephew Healthcare Ltd.
Healthcare House
101 Hessele Road
Hull, HU3 2BN
United Kingdom
Tel: 01482.222200
Fax: 01482.222211
npwtorders.uk@smith-nephew.com
24-hour Clinical Support Line:
UK 0845 120 0167

United States

Smith & Nephew, Inc.
Wound Management Division
970 Lake Carillon Dr.
St. Petersburg, FL 33716
Tel: 1 800-876-1261
Fax: 727-392-6914
customercare.largo@smith-nephew.com
1-866-998-NPWT (6798)

Made in USA for/Fabriqué aux Etats Unis Pour:
Smith & Nephew, Inc.,
970 Lake Carillon Drive, St. Petersburg, FL 33716

EC REP Smith & Nephew Medical Ltd.,
Hull HU3 2BN, England

°Trademark of Smith & Nephew
© Smith & Nephew 2012 www.smith-nephew.com